

examvue CR PRO

for Medical



Specifications:

- Cassette Sizes: 14x17",10x12",8x10"
**System Package includes two: 14x17" Cassette/IPs*
- System Set-up Options: Tabletop • Vertical wall mount
- Resolution: 100um/200um
- Data Capture:16-bits per pixel, 65,535 gray scale
- Dimensions (HxWxD): 5.3" x 17.5" x 31.3" (134 x 444 x 795 mm)
- Weight: 46lbs (21.5kg)

Workstation and Monitor

- Windows 7 Professional 32 Bit
- Intel Processor
- 500 GB HDD
- 4 GB RAM
- 16x CD/DVD +/- RW
- 21.5" LED Monitor

High
Resolution
Images

70
PPH

Acquire Software:

- DICOM Send,DICOM Print,CD/DVD Burning
- Multi Frequency Image Processing
- Zoom, Cropping, Mark, Window Level
- User Defined LUT
- Saving of measurements & annotations
- Embedded viewer

Viewing License:

- Distance Comparison
- Image Stitching
- Cobb Angle Measurements
- PACS connectivity for Easy Image Comparison
- Save Measurements and Annotations
- Prints on DICOM and Regular Printers
- DICOM Query/Retrieve /Send SCU/SCP
- DICOM Grayscale Print SCU
- View via Network

Regulatory approvals:

- FDA 510(k) Premarket Notification
- CE (0120)
- UL, cUL, FCC
- Health Canada Certification
- Not approved for mammography

Easy
to use
Software

Lightweight



*JPI Reserves the right to change specification without notice

JPI Healthcare Solutions Inc.
52 Newtown Plaza Plainview NY 11803
Toll Free: 888-574-2001
Tel: 516-513-1330
Fax: 516-513-1331
www.jpihealthcare.com

Authorized Dealer