



**Call Paul's X-Ray inc. (951)681-3241**

### Specifications:

- Cassette Sizes: 14x17", 10x12", 8x10"  
\*System Package includes two: 14x17" Cassette/IPs
- System Set-up Options: Tabletop • Vertical wall mount
- Resolution: 100um/200um
- Data Capture: 16-bits per pixel, 65,535 gray scale
- Dimensions (HxWxD): 5.3" x 17.5" x 31.3" (134 x 444 x 795 mm)
- Weight: 46lbs (21.5kg)

### Workstation and Monitor

- Windows 7 Professional 32 Bit
- Intel Processor
- 500 GB HDD
- 4 GB RAM
- 16x CD/DVD +/- RW
- 21.5" LED Monitor

High Resolution Images

70 PPH

### Acquire Software:

- DICOM Send, DICOM Print, CD/DVD Burning
- Multi Frequency Image Processing
- Zoom, Cropping, Mark, Window Level
- User Defined LUT
- Saving of measurements & annotations
- Embedded viewer

### Viewing License:

- Distance comparison
- Norberg angle and Cobb angle measurements
- Vertebral heart score (VHS)
- TPA Measurement (optional)
- PACS connectivity for easy image comparison
- Prints on both DICOM and off ice printers
- DICOM Send/Receive/Print
- Viewing via network/internet
- Image stitching

Easy to use Software

Lightweight

### Regulatory approvals:

- FDA 510(k) Premarket Notification
- CE (0120)
- UL, cUL, FCC
- Health Canada Certification

### Dental Upgrade:

- 10x14 Cassette housing and template
- Imaging Plate Starter Kit: Two(2) size 2 IP size  
One(1) x IP size 4

