

ImaSight 4600H DR

Chiropractic System

It's about *precision.*

With 16 million pixels and phenomenal HD images to present to your patients, razor-sharp images can help empower your practice, educate your patients and reflect the high standard of patient care you've set in the profession.

It's about *speed.*

From zero to diagnosis in less than 10 seconds: High-definition images from the ImaSight 4600H arrive in seconds, permitting immediate and faster diagnosis. And modifying the image is a snap - adjust contrast, brightness, and perspective on-the-fly and in real-time - *that's the digital revolution.*



It's about *your bottom line.*

With easy-to-use ImaSight Vision™ software annotation tools at your fingertips, you no longer have to explain what you're seeing to your patients, *you can show them* - giving them a strong reason to book their next appointment. Stop talking TO your patients - and start communicating WITH them.

It's about time... *to expect more.*

The reality is that film-based practices, just like the consumer film industry itself before it, is now giving way to the advantages only found in digital x-ray technology and advanced software.



Don't just make an image. **Make an impact™.**

ImaSightChiro.com | Toll Free: 888-707-4041 | info@imasightchiro.com



One
demonstration
can change
everything.

Knowledge is power: Make an informed decision by attending one of our weekly webinars or scheduling an online demonstration. Visit www.ImaSightChiro.com - and see for yourself the difference digital technology can make in your practice.



IMASIGHT

ImaSight 4600H DR Specs

Sensor

- Direct to Digital (DR) 16 Megapixels
- Image Acquisition time: 7 - 15 seconds
- Target Area: 14" x 17" / 35.56cm x 43.18cm
- Resolution: 4.6 lp/mm / Pixel Pitch: 108um
- DICOM compatible images
- ADC 16 bits
- Maintenance free

Acquisition Workstation PC

- Dell Inspiron 620 Intel® Core™ i3-2120 processor (3MB Cache, 3.30GHz)
- 4GB Memory
- Windows 7 Professional
- 1 TB Hard Drive
- External Backup Drive – 1 TB

High-Resolution Monitor

- Touchscreen
- 23" HD touch Screen
- 2.0 Megapixel [1980 x 1080] LCD
- Monitor Arm

Cables/Power Supply

- Ethernet cables (15'/4.57m + 3'/0.91m) + Ethernet isolator
- Ethernet Network Cable (14'/4.26m)
- 6 input surge suppressor power strip


ImaSight Vision Software

- One (1) Vision Capture Station License.
- Two (2) Vision Diagnostic Station Licenses.
- Includes a two (2) year subscription for all software updates and new releases, toll-free telephone and remote online software support.
- Additional ImASight Vision software licenses are available.

Warranty

- Five (5) years on sensor.
- Other components, such as computer, monitor, and related equipment carry manufacturer's warranties.

Show your patients powerful HD x-ray images in seconds - *anywhere* in your practice.



If a picture is worth a thousand words, imagine what an annotated x-ray is worth!

With ImASight in your practice, you can:

- Not only tell your patients your diagnosis - but show them. This leads to increased patient visits, increased profitability, and a fast ROI!
- Finally transform your x-ray room from a cost center to a profit center!
- Qualify for financing through our financial network - regardless of credit score!
- Adjust contrast, brightness, and perspective on-the-fly.
- Free yourself from maintenance worries - unlike our competitors' equipment - the 4600H does not contain moving parts.
- Benefit from a sensor that requires no special cooling and has no electronics exposed to x-rays.
- Say goodbye to film: Leave storage issues, messy and costly chemicals, and high doses of radiation where they belong - in the past.



Optional configuration. Large screen and laptop not included.

The ImASight4600H DR is part of the Infusion and Ultimate Systems - only available from ImASight. Learn more about these systems at www.ImASightChiro.com



"From day one of calling ImASight through its purchase and installation, my questions were handled thoroughly and with a personal touch. I researched all the other leading units and ImASight was simply the best." - Dr. Roger Borbon