



ImaSight 4600 Veterinary Solutions

IMASIGHT 4600 – TARGET AREA - 14 X 17

- Resolution 4.6 lp/mm • Pixel Pitch:108 µm
- ADC 16 bits

ACQUISITION WORKSTATION: DELL OPTI 360 CORE DUO, 2.66 GHZ

- 2 GB RAM
- 320 GB Hard Drive
- External Backup Drive – 500 GB

HIGH-RESOLUTION MONITOR

- Touch screen
- 20" viewable area
- 2.0 MP, 1600 x 1200 LCD

POWER SUPPLY (WITH MEDICAL POWER CORD)

- 6 input surge suppressor power strip

CABLES

- Imasight 4600 Sensor 10/100 BASE Ethernet cables + Ethernet isolator
- DB-9 serial cable (15')
- Network Ethernet Cable (14')

IMASIGHT VISION SOFTWARE

- Product package includes an additional license for Imasight Vision Software for use on one clinic diagnostic workstation.
- Additional Imasight Vision software licenses are available.
- Two (2) year subscription for all software updates and new releases, which includes toll-free telephone and remote online software support at no extra charge. Extension of this subscription is available at \$1,000 per year (price subject to change). Software support is available on an hourly rate basis without subscription, if required.

WARRANTIES

- Five (5) years on sensor
- Two (2) years on parts
- Other components, such as computer, monitor, and related equipment carry manufacturers' warranties.



Technical Support Hotline
1.866.989.7979
Weekdays 8am-8pm EST
Saturday 9am-3pm EST

**YOUR PRACTICE.
OUR FOCUS.**

925 Carriere Blvd, Suite 10
Gatineau QC J8Y 6W5 CANADA

PHONE: 819.777.7979
FAX: 613.822.5195

IMASIGHT.COM

Distributed by



190 Boul. St-Elzéar O., Laval (Qc) H7L 3N3

Tel : 450 669-8985 □ 1 800 361-0877

Fax : 450 669-9532 □ 1 888 935-7001

www.medicaltronik.ca



ImaSight 4600 Veterinary Table Specifications

High Frequency Generator

- 30 kW
- 125 kVp
- Membrane Console
- 400 mA
- Anatomical Programming
- Power Supply Cabinet: 27" x 14" x 45"

Table

- Fixed 40 inch SID Tube Stand
- High Voltage Cables (15')
- 5 foot Table Top
- Exposure Foot Pedal
- Manual Collimator
- Four-way Floating Top Table

X-Ray Tube

- 0.6/1.5 Focal Spots
- 200,000 Hu

Warranties

- X-Ray System: 5 years
- X-Ray Tube: 1 year





IMASIGHT

ImaSight Solution Offering

ImaSight Support Team

- ImaSight has a team of qualified and trained people who are taking care of our customers directly. Hotline number: (866) 989-7979 Email address: support@imasight.com
- ImaSight remotely supports installers during each installation. Quality Assurance tests are performed at every site. All units are tested and pre-calibrated before shipping.
- A site survey is sent to the clinics in order to plan for installations and eliminate mistakes. Staff Training is planned early to reduce downtime.
- ImaSight offers webinars to dealers on our easy-to-learn software.

Updates and Upgrades

- ImaSight is constantly improving its software with its own team of programmers. Suggestions from dealers are welcome and are integrated within our core software development.
- Our GUI software evolves with the market while keeping the core very reliable. Software updates are released every six months.

Retrofits available

- ImaSight offers, at no charge, retrofits for a wide range of manufacturer's tables. Call to inquire about our retrofit kits.

Detector Only Pricing

- ImaSight can sell the detector alone to integrators who prefer to integrate the sensor within the table of their choice.
- Call for pricing : (866) 217-8183 or sales@imasight.com



Distributed by



190 Boul. St-Elzéar O., Laval (Qc) H7L 3N3

Tel : **450 669-8985** □ **1 800 361-0877**
Fax : **450 669-9532** □ **1 888 935-7001**

www.medicaltronik.ca

**YOUR PRACTICE.
OUR FOCUS.**

925 Carrière Blvd, Suite 10
Gatineau QC J8Y 6W5 CANADA

PHONE: 819.777.7979
FAX: 613.822.5195

IMASIGHT.COM





IMASIGHT

ImaSight Network Solution

ImaSight Vision
Mobile Viewing Station #1



ImaSight Vision
Viewing Station #1



Your Clinic Network



ImaSight Vision
Viewing Station #2



X-Ray Room
ImaSight
Capture Workstation



...

You can add as many
ImaSight Viewing Stations
as you want.

We also offer PACS network solutions.

Distributed by



190 Boul. St-Elzéar O., Laval (Qc) H7L 3N3

Tel : 450 669-8985 □ 1 800 361-0877

Fax : 450 669-9532 □ 1 888 935-7001

www.medicaltronik.ca 

YOUR PRACTICE.
OUR FOCUS.

925 Carrière Blvd, Suite 10
Gatineau QC J8Y 6W5 CANADA

PHONE: 819.777.7979
FAX: 613.822.5195

IMASIGHT.COM



IMASIGHT

PACS Network Solution

3rd Party Software
Viewing Station #1



Your Clinic Network



3rd Party Software
Viewing Station #2




3rd Party Software
Viewing Station #3



X-Ray Room
**ImaSight
Capture Workstation**



PACSDICOM Server



Hybrid Solutions - ImaSight Vision Workstations and 3rd Party Software Workstations on the same network - are also possible.

Distributed by



190 Boul. St-Elzéar O., Laval (Qc) H7L 3N3

Tel : 450 669-8985 □ 1 800 361-0877
Fax : 450 669-9532 □ 1 888 935-7001

www.medicaltronik.ca



**YOUR PRACTICE.
OUR FOCUS.**

925 Carrière Blvd, Suite 10
Gatineau QC J8Y 6W5 CANADA

PHONE: 819.777.7979
FAX: 613.822.5195

IMASIGHT.COM